



Multifunctional Ultrasonic Cleaner

LMFU-A12

Multifunctional Ultrasonic Cleaner LMFU-A12

Multifunctional Ultrasonic Cleaner LMFU-A12 is a single 15 L multifunctional tank with heat and ultrasonic action providing deep cleansing, removing contamination even from blind holes and internal surfaces. Its multifunctionality includes power adjusting mode, pulse or degas control and variable ultrasonic intensity control. Equipped with degassing and degreasing function for removal of dissolved gasses and cleansing of the tank.

Features

- ❖ 15 L capacity tank
 - ❖ Integrated digital circuit for timer and heater
 - ❖ Temperature range : 0 °C ~ 100 °C
 - ❖ Power adjustable mode
 - ❖ Cleaning duration : 1 ~ 99 mins (can be kept 'ON' always)
 - ❖ 4 working modes available (Power , Pulse , Turbo , Degas)
 - ❖ One key access for auto stand-by , sleep and wake up mode
 - ❖ On power cut, auto memory save option
 - ❖ Efficient cleaning by distilled water
 - ❖ Equipped with degassing and degreasing function
 - ❖ Precise digital processed temperature controls
 - ❖ Optimum durability under prolonged and heavy usage
-

Application

Ideal for cleaning irregular and complex instruments, lenses and other optical parts, dental and surgical instruments, electronics. It is also used in extraction and to degas certain analyzed stuff in laboratories

Multifunctional Ultrasonic Cleaner LMFU-A10

Specifications

Model No.	LMFU - A12
Capacity	15 L
Tank size (L x W x H)	330 x 300 x 150 mm
Overall size (L x W x H)	360 x 330 x 357 mm
Ultrasonic power	0 - 360 W
Heater power	400 W
Temperature set	0 °C ~ 100 °C
Timer	1 ~ 99 mins (optional : always ON)
Construction material	Stainless Steel
Drainage	Yes
Voltage	AC 220 V / 50 Hz