

Mechanical Heating Ultrasonic Cleaner

LMUC-A18



Mechanical Heating Ultrasonic Cleaner LMUC-A18 delivers high efficiency functionality and features. This compact model with 30 L tank capacity and 600 W ultrasonic power provides deep cleansing action removing contamination even from blind holes and internal surfaces. Its key feature is the variance of using it at 28 KHz for gross rinse to 40 KHz for intensive rinse Equipped with degassing and gettering function for removal of dissolved gasses and efficient cleansing of oil and wax

Features

- ◆ 30 L capacity tank
- ◆ Integrated mechanical circuit for timer and heater
- ◆ Temperature range : 20 °C ~ 80 °C (mechanical)
- ◆ Cleaning duration : 0 ~ 20 mins
- ◆ Efficient cleaning can be done by distilled water
- ◆ Optimum durability under prolonged and heavy usage
- ◆ Equipped with degassing and gettering function

Application

Ideal for cleaning irregular and complex instruments, lenses and other optical parts, dental and surgical instruments, electronics. It is also used in extraction and to degas certain analyzed stuff in laboratories

Specifications

Model No.	LMUC-A18
Capacity	30 L
Tank size (L x W x H)	500 mm x 300 mm x 200 mm
Machine size (L x W x H)	550 mm x 330 mm x 380 mm
Ultrasonic power	600 W
Heater power	500 W
Temperature set	20 °C ~ 80 °C
Digital timer (mechanical)	0 ~ 20 mins
Control mode	Mechanical
Construction material	Stainless Steel
Voltage	AC 220 V / 50 Hz