



POWER ADJUSTABLE ULTRASONIC CLEANER

LPAU-A10

Power Adjustable Ultrasonic Cleaner LPAU-A10

Power Adjustable Ultrasonic Cleaner LPAU-A10 is the most reliable instrument with 6.5 L of tank capacity and cleaning duration of 1 - 99 mins for efficient intensive cleaning of instruments and products making them contamination, rust and germ free. Its key feature is the variance of using it at 28 KHz for gross rinse to 40 KHz for intensive rinse. The power can be adjusted which provides the user with a precise control from zero watts to maximum power.

Features

- ◆ 6.5 L capacity tank
 - ◆ Temperature range : 0 °C ~ 80 °C
 - ◆ Cleaning duration : 1- 99 mins
 - ◆ 28 KHz for gross rinse to 40 KHz for intensive rinse
 - ◆ One key access for auto stand-by , sleep and wake up mode
 - ◆ Integrated digital circuit for timer and heater
 - ◆ Precise digital processed temperature controls
 - ◆ Efficient cleaning by distilled water
 - ◆ Optimum durability under prolonged and heavy usage
 - ◆ Equipped with degassing and degreasing function
-

Application

Ideal for cleaning irregular and complex instruments, lenses and other optical parts, dental and surgical instruments, electronics. It is also used in extraction and to degas certain analyzed stuff in laboratories

Specifications

Model No.	LPAU - A10
Capacity	6.5 L
Tank size (L x W x H)	300 x 155 x 150 mm
Overall size (L x W x H)	330 x 180 x 337 mm
Ultrasonic power	0 ~ 180 W
Heater power	300 W
Temperature set	0 °C ~ 80 °C
Timer	1-99 mins
Construction material	Stainless Steel
Drainage	Yes
Voltage	AC 220 V / 50 Hz