



# ANAEROBIC JAR

## LAJ-B20

Anaerobic Jar LAJ-B20 is a hermetic type jar comprising transparent jar and lid with 25 R petri-dish holder. Features PMMA material with extremely smooth surface and excellent air tight performance. Designed with strong jar-lid clamp construction and O-ring sealing, this anaerobic jar is adopted with quick snap-shut coupling, enables one swift move connection to system. With colored stainless-steel clamp and power electro-static painting, it is shock and impact resistant with durability and convenience for experimental management.

### Features

- ❑ Hermetic type jar with transparent jar and lid
- ❑ 25 R petri-dish holder with single stack of 12 dishes ( $\phi 90$  to 100 mm)
- ❑ PMMA material with extremely smooth surface and excellent air tight performance
- ❑ Strong jar-lid clamp construction and O-ring sealing, for outstanding performance
- ❑ Adopted with quick snap-shut coupling, enables one swift move connection to system
- ❑ Colored stainless-steel clamp and power electro-static painting
- ❑ Durability and convenience for experimental management
- ❑ Shock and impact resistant, easy to clean
- ❑ Lightweight design with convenient portability

### Application

Anaerobic Jars are used for anaerobic incubation (anaerobiosis) of bacterial cultures for research and testing across the fields of microbiology, biotechnology, pathology, virology and CDC.

## Specifications

| Model                                      | LAJ-B20  |
|--|--|
| Type                                       | Hermetic type                                    |
| Time taken (to attain anaerobic condition) | 4 hours  |
| O-ring                                     | Yes  |
| Contents                                   | Single stack of 12 dishes ( $\phi$ 90 to 100 mm) |
| Petri dish holder                          | 25 R   |
| Material                                   | PMMA material                                    |
| Clamp                                      | Colored stainless-steel clamp                    |
| Quick snap-shut coupling                   | Yes  |
| Packaging dimensions (L×W×H)               | 245×215×395 mm                                   |
| Gross weight                               | 2.5 kg   |