

Manual Dual Door

Toxic Cabinet

LTXC-A10



170 L Toxic Cabinet LTXC-A10

170 L Toxic Cabinet LTXC-A10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 170 L
 - Dual walled manual
 - 38 mm insulating air space
 - 180° fully opening door
 - Three point latch and manual lock
 - Two adjustable shelves
 - Adjustable spill retaining shelves
-

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications

Model No	LTXC-A10
Capacity	170 L
Door type	Dual door manual
Adjustable shelves	2
Shelves loading capacity	100 kgs
Dimension (H x W x D)	1650 x 1090 x 460 mm
Weight	148 kgs