

Manual Dual Door
Corrosive Cabinet
LCOC-A1 SERIES



114 L Corrosive Cabinet LCOC-A10

114 L Corrosive Cabinet LCOC-A10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 114 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

170 L Corrosive Cabinet LCOC-A11

170 L Corrosive Cabinet LCOC-A11 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 170 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications

Model No	LCOC-A10	LCOC-A11
Capacity	114 L	170 L
Door type	Dual door manual	Dual door manual
Adjustable shelves	1	2
Shelves loading capacity	50 kgs	100 kgs
Dimension (H x W x D)	1120 x 1090 x 460 mm	1650 x 1090 x 460 mm
Weight	109 kgs	151 kgs

227 L Corrosive Cabinet LCOC-A12

227 L Corrosive Cabinet LCOC-A12 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 227 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

340 L Corrosive Cabinet LCOC-A13

340 L Corrosive Cabinet LCOC-A13 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 340 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications

Model No	LCOC-A12	LCOC-A13
Capacity	227 L	340 L
Door type	Dual door manual	Dual door manual
Adjustable shelves	2	2
Shelves loading capacity	100 kgs	100 kgs
Dimension (H x W x D)	1650 x 860 x 460 mm	1650 x 1090 x 860 mm
Weight	151 kgs	186 kgs