

**Manual Single Door**  
**Corrosive Cabinet**  
LCOC-B1 SERIES



### 15 L Corrosive Cabinet LCOC-B10

15 L Corrosive Cabinet LCOC-B10 is a dual walled manual combustion cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- Capacity - 15 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

### 45 L Combustible Cabinet LCOC-B11

45 L Combustible Cabinet LCOC-B11 is a dual walled manual combustion cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- Capacity - 45 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

#### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specification

Model No	LCOC-B10	LCOC-B11
Capacity	15 L	45 L
Door type	Single door manual	Single door manual
Adjustable shelves	1	1
Shelves loading capacity	50 kgs	50 kgs
Dimension ( H x W x D )	560 x 430 x 430 mm	890 x 590 x 460 mm
Weight	50 kgs	60 kgs

## 83 L Combustible Cabinet LCOC-B12

83 L Combustible Cabinet LCOC-B12 is a dual walled manual combustion cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features three adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

## Features

- Capacity - 45 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Three adjustable shelves
- Adjustable spill retaining shelves

## Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

### Specification

Model No	LCOC-B12
Capacity	83 L
Door type	Single door manual
Adjustable shelf	3
Shelves loading capacity	50 kgs
Dimension ( H x W x D )	1650 x 600 x 460 mm
Weight	101 kgs