



PARAFFIN WAX DISPENSER LPWD-A10

Paraffin Wax Dispenser LPWD-A10 is a compact and automated unit adopted with computer control and intelligent digital temperature control, offers fast and reliable heating. With advanced heating element and temperature sensors, features 10000 ml capacity with 0 to 100°C of continuously adjustable temperature range. Designed with Iron coated housing, stainless steel interior, LED display and memory and automatic restoration function, it can automatically save last operation setting and allows automatic saving of preset temperature after operation.

Features

- ❑ A compact and automated unit with computer control, easy and convenient operation
- ❑ Intelligent digital temperature control, precise, accurate and reliable heating
- ❑ Advanced heating element and temperature sensors, fast and uniform heating
- ❑ 10000 ml capacity with 0 to 100°C of continuously adjustable temperature range
- ❑ Iron coated housing with stainless steel interior, easy maintenance and long lifespan
- ❑ LED display, allow convenient setting and monitoring of parameters
- ❑ Automatically saves the last operational setting, reproducible results
- ❑ Memory and automatic restoration function, allows automatic saving of preset temperature after operation
- ❑ High-efficient and stable unit for having molten paraffin available at all times

Application

Paraffin Wax Dispenser is used to heat paraffin wax to a liquid state and provides convenient wax on tap across histology laboratories, forensics, medicine, anatomy etc.

Specifications

Model	LPWD-A10
Capacity	10000 ml
Temperature range	0 to 100°C (continuously adjustable)
Temperature control precision	± 1°C
Temperature control	Intelligent digital temperature control
Display	LED
Power supply	AC 220/110 V ± 10%, 50/60 Hz
Power	1200 W
Dimension (W×D×H)	355×410×540 mm
Net weight	17 kg