



PHARMACY REFRIGERATOR

LPRF-A15

Pharmacy Refrigerator LPRF-A16 is a digital microprocessor controlled temperature refrigerator with holding capacity of 360 L, 2 to 8°C operating temperature range, and adjustable storage shelves. Features Inbuilt full automatic defrost system with CFC-Free R600a refrigerant. The electrical heater decreases condensation on glass door in high humidity environment. With user-friendly operation mode, refrigerator is designed with audible and visual alarm system for different conditions.

Features :

- ❑ Accurate temperature control system maintains inner temperature in 2 to 8 °C
- ❑ High-tech integrated sensors to display and control temperature
- ❑ Auto-defrost to remove moisture on cooling surface
- ❑ Removable outdoor double seal design: better temperature maintain and energy saving
- ❑ Microprocessor controlled forced-air cooling system
- ❑ Digital temperature displays with 0.1°C resolution
- ❑ CFC-Free R600a refrigerant and Shelves adjustable design
- ❑ Epoxy coated exterior anti-bacteria and easy cleaning
- ❑ Caster with break design for easy moving and lock
- ❑ Interior LED light for easy observation
- ❑ Electrical heater to decrease condensation on door in high humidity environment
- ❑ Different Alarm conditions High/Low Temperature; Power Failure; Sensor Error; Door Ajar, Low Voltage in Backup Battery

Application :

The pharmacy refrigerator finds application across medical research centers, drug stores, pharmacy companies, hospitals, epidemic prevention centre's to store medicines, testing kits, reagents, vaccines etc.

Specifications :

Model	LPRF-A15
Cabinet Type	Upright
Gross Capacity	360 L
Climate Class	10 to 32°C
Cooling Type	Forced Air Cooling
Defrost Mode	Automatic
Refrigerant	CFC-Free R600a
Door	Glass Door
Noise	50 dB
Temperature Range	2 to 8°C
Uniformity	± 0.1°C
Temperature Control	Microprocessor
Display	LED digital display
Power Supply	220 V, 50 Hz and 138 W
Current	0.88 A
Power Consumption	1.78 kWh/24h
Internal Dimension (W×D×H)	540 × 425 × 1380 mm
External Dimension (W×D×H)	640 × 582 × 1880 mm
Packing Dimension (W×D×H)	710 × 655 × 1960 mm
Net Weight	92 kg
Gross Weight	100 kg