



Platelet Shaker

LPS-A1 Series

Platelet Shaker LPS-A1 Series is made of high-grade steel with stainless steel trays for easy cleaning and durability. Supplied with multiple trays to accommodate 12 to 48 platelet bags, offers low maintenance with quiet operation. Platelet shakers specifically designed to meet the demands of various clinical laboratories.

Features :

- ❑ Platelet bags can be loaded easily on the stainless steel basket
- ❑ Mesh Design Provides air circulation
- ❑ Exhibits a steady and smooth 360° movement
- ❑ Attached with a tachometer for easy monitoring of speed
- ❑ Robust and practical design makes it economical choice

Application :

Platelet shakers are widely used in blood banks and clinical research laboratories.

Specifications:

Model	LPS-A10	LPS-A11
Shaking Speed Range	0 to 9 rpm	0 to 9 rpm
Capacity	12 bags	24 bags
Overall Dimension	370 x 255 x 390 mm	510 x 255 x 290 mm
Power Supply	220 V, 60 Hz	220 V, 60 Hz

Model	LPS-A12	LPS-A13
Shaking Speed Range	0 to 9 rpm	0 to 9 rpm
Capacity	36 bags	48 bags
Overall Dimension	650 x 255 x 290 mm	790 x 255 x 290 mm
Power Supply	220 V, 60 Hz	220 V, 60 Hz