

Portable Coffee Refractometer

LCR-A10



Portable Coffee Refractometer LCR-A10

Portable coffee refractometer LCR-A10 is a microprocessor controlled pocket sized optical instrument that provides instant measurement readings. LCD screen displays measurement data in one view. Automatically turns off in 1 minute of non-use prolonging it's battery life. It is designed for measuring brix and TDS of coffee.

Features

- Digital meter displays easy to read results
- Standardized Brix scale for easy read and accuracy
- Automatic and zero calibration with distilled or tap water
- Abrasion and corrosion-resistant coating on sapphire prism
- Auto power off (after 1 minute)

Applications

Used in laboratories, research institutes, beverage industry, food industry, etc.

Specifications

Model	LCR-A10
Brix range	0.00 ~ 30 %
TDS range	0.00 ~ 25.00
Temperature	10 °C ~ 100 °C
Battery	1.5 V
Dimensions	121 × 58 × 25 mm
Weight	90 grams