



Rapid Low Temperature Circulating Bath

LRTB-A20

www.labtron.com
info@labtron.com

Rapid Low Temperature Circulating Bath LRTB-A20

Rapid Low Temperature Circulating Bath LRTB-A20 is a microcomputer controlled 5 L bath with refrigeration system and a temperature range of - 40 °C to 200 °C. Equipped with dual red and green LED windows with upper window measuring red value and lower measuring green value. Features internal and external cycles with external cycles discharging thermo-static liquid outside the bath raising the temperature of the external vessel.

Features

- ▶ Volume - 5 L
- ▶ Low temperature cooling pump
- ▶ Microcomputer controlled temperature
- ▶ Temperature and time LCD display
- ▶ Built-in liquid level monitor
- ▶ Refrigeration unit with relay , protector and capacitor
- ▶ Anti-corrosive circulatory system

Specifications

Model No	LRTB-A20
Volume	5 L
Temperature range	- 40 °C to 200 °C
Refrigerating capacity	1322 W
Pump flow	20 L / min
Temperature stability	± 0.1 °C
Head of delivery	10 m
Temperature control tier number	32
Heating power	800