



# Rapid Low Temperature Circulating Bath

LRTB-A73

## Rapid Low Temperature Circulating Bath LRTB-A73

Rapid Low Temperature Circulating Bath LRTB-A73 is a microcomputer controlled 30 L bath with cold ring liquid circulation and cryogenic liquid with a temperature range of -40 °C. Refrigeration unit with low temperature cooling pump facilitates consistent high pressure and maintenance of low temperature under working conditions of reaction experiments. Equipped with anti-corrosive circulatory system preventing traffic flow, high pressure and low pump body temperature.

### Features

- ▶ Volume - 30 L
- ▶ Low temperature cooling pump
- ▶ Microcomputer controlled temperature
- ▶ Temperature and time LCD display
- ▶ Built-in liquid level monitor
- ▶ Refrigeration unit with relay, protector and capacitor
- ▶ Anti-corrosive circulatory system

### Specifications

Model No	LRTB-A73
Volume	30 L
Temperature range	-40 °C
Refrigerating capacity	11808 W
Pump flow	20 L / min
Temperature stability	±0.5 °C
Head of delivery	10 m