

Self Close Dual Door

Flammable Storage Cabinet

LFSC-C1 Series



170 L Flammable Storage Cabinet LFSC-C10

170 L Flammable Storage Cabinet LFSC-C10 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 170 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications

Model No	LFSC-C10
Capacity	170 L
Door type	Dual door self-closing
Adjustable shelves	2
Shelves loading capacity	100 kgs
Dimension (H x W x D)	1650 x 1090 x 460 mm
Weight	148 kgs

227 L Flammable Storage Cabinet LFSC-C11

227 L Flammable Storage Cabinet LFSC-C11 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 227 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications

Model No	LFSC-C11
Capacity	227 L
Door type	Dual door self-closing
Adjustable shelves	2
Shelves loading capacity	100 kgs
Dimension (H x W x D)	1650 x 860 x 860 mm
Weight	173 kgs

340 L Flammable Storage Cabinet LFSC-C12

340 L Flammable Storage Cabinet LFSC-C12 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features

- Capacity - 340 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications

Model No	LFSC-C12
Capacity	340 L
Door type	Dual door self-closing
Adjustable shelves	2
Shelves loading capacity	100 kgs
Dimension (H x W x D)	1650 x 11090 x 860 mm
Weight	186 kgs