

**Single Door Self-Close**

# Flammable Storage Cabinet

**LFSC-D1 Series**



### 15 L Flammable Storage Cabinet LFSC-D10

15 L Flammable Storage Cabinet LFSC-D10 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Door is reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelf for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- Capacity - 15 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelf

#### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

#### Specifications

Model No	LFSC-D10
Capacity	15 L
Door type	Single door self-closing
Adjustable shelf	1
Shelves loading capacity	50 kgs
Dimension ( H x W x D )	560 x 430 x 430 mm
Weight	30.5 kgs

### 45 L Flammable Storage Cabinet LFSC-D11

45 L Flammable Storage Cabinet LFSC-D11 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Door is reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelf for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- Capacity - 45 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelf

#### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

#### Specifications

Model No	LFSC-D11
Capacity	45 L
Door type	Single door self-closing
Adjustable shelf	1
Shelves loading capacity	50 kgs
Dimension ( H x W x D )	890 x 590 x 460 mm
Weight	59 kgs