



5 L DISTILLER



10 L DISTILLER



20 L DISTILLER

LSWD-A2 Series

Stainless Steel Water Distiller

Stainless Steel Water Distiller LSWD-A2 series

Stainless Steel Water distiller LSWD-A2 series are single reservoir distillation units with automatic flow-control. Stainless steel construction aids to produce germ-free distillate. Narrow piping connects heater column to heat exchanger for optimal distillate quality. Light switch indicator signifies the operating status of the unit. Branched valve control guides flow passage for safe use.

Features

- USP grade pyrogen-free distillate
- Shunt valve regulates inlet and outlet flow
- Gas/stench exhaustion vent
- Corrosion-proof reservoir
- Auto-cutoff water supply for exceeded water levels
- Control valve for inlet water flow

Applications

Used in laboratory, organic chemistry lab, fermentation and medical industry to produce industrial grade highly treated and distilled water suitable for laboratory usage.

Specification

MODEL NO	LSWD - A20	LSWD - A21	LSWD - A 22
Capacity	5 L/hr	10 L/hr	20 L/hr
Water consumption	1:10	1:10	1:10
Configuration	Automated control	Automated control	Automated control
Water type	Regular water	Regular water	Regular water
Power	220 V ; 2.5 kW	380 V ; 2.5 Kw	380 V ; 4.5 Kw
	Single phase	Triple phase	Triple phase
Dimensions	365 × 368 × 845 mm	365 × 368 × 845 mm	365 × 368 × 945 mm
Weight	9 kgs	9 kgs	10 kgs