



# TABLE TOP AUTOCLAVE

## LT TA-A20

## Table top autoclave LTTA-A20

Table top autoclave LTTA-A20 is a robust and compact, space efficient unit with 20 L of capacity and a sterilizing temperature of 134 °C. It is equipped with a key feature of 4 - 6 mins rapid sterilization along with safety features like automatic shutoff preventing overheating of chamber, automatic cold air discharge and auto protection against raised temperature and pressure. Integrated with steam-water inner circulation system with steam discharge making the environment for sterilization clean and dry

### Features

- Digital LCD display to monitor working status and parameters
- Working temperature : 134 °C
- 4 - 6 mins rapid sterilization
- Automated drying function
- Equipped with interior circulate system
- Inner steam circulatory system
- Automatic cold air discharge from system chamber ensures dry and clean environment
- Door lock feature : Door cannot be opened till the pressure inside the chamber doesn't fall to 0.027 MPa
- Self-inflating leak proof chamber
- Auto exhaustion of cool air
- Safe door lock system ensuring safety / protection during sterilization
- Automatic power cut off with alarm with low water levels indication
- Stainless steel sterilizing baskets

### Application

Used to sterilize culture media, instruments, dressings, intravenous equipments, rubber closure of glass equipments used in aseptic techniques and numerous objects that can withstand high temperature and pressure

## Specification

Model	LTTA-A20
Capacity	24 L
Number of sterilizing trays	3
Sterilization Temperature	134 °C
Sterilization Pressure	0.22 Mpa
Timer Range	0-60 Minutes
Heat Average	$\leq \pm 1$ °C
Chamber Dimension ( $\Phi \times D$ )	270×410 mm
Net Weight	18 Kgs
Gross Weight	20 Kgs
Power Supply	AC 220 V 50 Hz/1.5 KW